verifi[®]
ORVOSOKNAK

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1. American College of Obstetricians & Gynecologists. ACOG Practice Bulletin 77. Screening for Fetal Chromosomal Abnormalities. *Obstet Gynecol.* 2007;109:217–228.
2. American College of Obstetricians & Gynecologists. Committee Opinion No. 545. Noninvasive Prenatal Testing for Fetal Aneuploidy. *Obstet Gynecol.* 2012;120:1532–1534.

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
Electronically signed and dated on 12/6/2013 5:08:12 PM
William K. Seltzer, PhD, FACMG

PROVIDER INFORMATION
Elizabeth Blackwell
Hortelown Hospital
21 Elm Street
New York, NY 10007

Report Copy To: J. Smith
Phone: (123) 456-7899

| TEST | PRELIMINARY RESULT |
|---------------|------------------------|
| Chromosome 21 | No aneuploidy detected |
| Chromosome 18 | No aneuploidy detected |
| Chromosome 13 | No aneuploidy detected |

Comments: Clinical correlation is desired, chorionic villus sampling

Order ID: 308177

Test Method: Nuclear acid extraction, DNA sequencing using next generation sequencing technology, and analysis of sequencing results to determine test aneuploidy.
Limitations of Test: This test is designed to detect chromosome aneuploidies of chromosomes 21, 18, and 13, X and Y. The test is not validated for use in pregnancies with gestational age of at least 10 weeks. These results do not eliminate the possibility that this pregnancy is associated with other chromosomal abnormalities, both detected and not detected. A negative test result does not provide the presence of a normal pregnancy. The presence or absence of aneuploidy cannot be determined. The presence or absence of aneuploidy cannot be determined. The presence or absence of aneuploidy cannot be determined. There is a small possibility that the test results might not reflect the chromosomal status of the fetus. The results of this test should be used in conjunction with other clinical information. See the limitations of this test on the Verifi prenatal test performance data.

Prenatal Aneuploidy Test Report

REPORT DATE AND TIME
04.Dec.2013 | 5:00 PM

Electronically signed and dated on 12/6/2013 4:58:58 PM
William K. Seltzer, PhD, FACMG

PROVIDER INFORMATION
Elizabeth Blackwell
Hortelown Hospital
21 Elm Street
New York, NY 10007

Report Copy To: J. Smith
Phone: (123) 456-7899 | Fax: (123) 456-7789

Prenatal Aneuploidy Test Report

ANEUPLOIDY DETECTED

See Below

REPORT DATE AND TIME
04.Dec.2013 | 4:58 PM

PROVIDER INFORMATION

Elizabeth Blackwell
Hortelown Hospital
21 Elm Street
New York, NY 10007

Report Copy To: J. Smith
Phone: (123) 456-7899 | Fax: (123) 456-7789

PATIENT INFORMATION

Name: Susan O Brown
DOB: 01 Jul 1990 | Sex: Female
Medical Record/Patient ID: 12356
Gestational Age at Draw: 10 weeks | 10
Indication: Positive Serum Screen
Genetic Counselor: Jennifer Jones
Client Sample ID:

PRENATAL ANEUPLOIDY TEST RESULTS – Singleton Pregnancy

| TEST | RESULT | INTERPRETATION |
|-----------------|------------------------|---|
| Chromosome 21 | ANEUPLOIDY DETECTED | Results consistent with trisomy for chromosome 21 |
| Chromosome 18 | No aneuploidy detected | Results consistent with diploid chromosome 18 |
| Chromosome 13 | No aneuploidy detected | Results consistent with diploid chromosome 13 |
| Sex Chromosomes | No aneuploidy detected | Results consistent with two sex chromosomes (XY) |

Comments: Genetic counseling is recommended. Clinical correlation with ultrasound findings, history, and other screening tests is indicated. If definitive diagnosis is desired, chorionic villus sampling or amniocentesis is necessary.

SAMPLE ID

| | | | |
|------------------|---------------------|------------------------|---------------------------|
| Order ID: 308177 | Sample ID: 308177-1 | Draw Date: 29.Nov.2013 | Receipt Date: 29.Nov.2013 |
|------------------|---------------------|------------------------|---------------------------|

verifi prenatal test performance data*

| Chromosome | N | Frequency | %N | Frequency | %N |
|------------|-----|-------------|------|-------------|------|
| 21 | 500 | 149 (29.8%) | 14.9 | 100 (20.0%) | 10.0 |
| 18 | 500 | 87 (17.4%) | 8.7 | 100 (20.0%) | 10.0 |
| 13 | 500 | 87 (17.4%) | 8.7 | 100 (20.0%) | 10.0 |
| Monosomy X | 500 | 10 (2.0%) | 2.0 | 100 (20.0%) | 10.0 |
| XXY | 500 | 10 (2.0%) | 2.0 | 100 (20.0%) | 10.0 |
| XYY | 500 | 10 (2.0%) | 2.0 | 100 (20.0%) | 10.0 |

*See the limitations of this test on the Verifi prenatal test performance data.

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CLIA

Clinical Laboratory Improvement Amendments

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Illumina Inc.